FLEX
Monoclonal Mouse
Anti-Human
Progesterone Receptor
Clone PgR 636
Ready-to-Use
(Dako Autostainer/Autostainer Plus)

English

Code IS068

Intended use	For in vitro diagnostic use. FLEX Monoclonal Mouse Anti-Human Progesterone Receptor, Clone PgR 636, Ready-to-Use, (Dako Autostainer/Autostainer Plus), is intended for use in immunohistochemistry together with EnVision™ FLEX+, High pH visualization kit together with Autostainer/Autostainer Plus instruments to semi-quantitatively detect human progesterone receptor in formalin-fixed, paraffin-embedded human breast carcinoma. This antibody labels progesterone receptor-positive cells and is useful in the assessment of progesterone receptor status in human breast carcinomas. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.
Synonym for antigen	PR
Summary and explanation	The role of steroid hormone receptors in breast cancer is well-known (1,2). The absence of estrogen receptor ER and progesterone receptor PR predicts early recurrence and poor survival of breast cancer patients (3-6). Also, the presence of ER and PR in tumors predicts the potential for benefit from tamoxifen and other endocrine-related therapies. Measurement of ER and PR can be determined semi-quantitatively using IHC or quantitatively using DCC or EIA. Correlation between the semi-quantitative and quantitative evaluations of PR have ranged from 73 to 91% depending on the laboratory and antibody used (7-9). Refer to Dako's <i>General Instructions for Immunohistochemical Staining</i> or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.
Reagent provided	Ready-to-use monoclonal mouse antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L NaN₃. <u>Clone</u>: PgR 636 (10) <u>Isotype</u>: IgG1, kappa
Immunogen	Formalin-fixed recombinant full-length A-form of human progesterone receptor (10)
Specificity	Anti-PR, PgR 636 has been demonstrated to react with the PR-A and PR-B forms by Western blot of whole cell extracts and reacts with both free and hormone-bound PR (10). The epitope has been mapped to the amino terminal domain shared by PR-A and PR-B.
Precautions	- For professional users. - This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing. - As with any product derived from biological sources, proper handling procedures should be used. - Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin. - Unused solution should be disposed of according to local, State and Federal regulations.

Storage	Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.
Specimen preparation including materials required but not supplied	<u>Paraffin sections:</u> The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm. <u>Pre-treatment:</u> Pre-treatment of formalin-fixed, paraffin-embedded tissue sections is required. Optimal results are obtained by pretreating tissues with HIER using diluted EnVision™ FLEX Target Retrieval Solution, High pH (50x) (Codes K8012/K8004). Deparaffinization, rehydration and epitope retrieval can be performed in Dako PT Link (Code PT100/PT101). For details, please refer to PT Link User Guide. The following parameters should be used for PT Link: Pre-heat temperature: 65 °C; epitope retrieval temperature and time: 97 °C for 20 (±1) minutes; cool down to 65 °C. Remove slide rack from PT tank and immediately dip slides in jar/tank (e.g., PT Link Rinse Station (Code PT109)) containing diluted room temperature EnVision™ FLEX Wash Buffer (20x) (Code K8007). Leave slides in Wash Buffer for 1-5 minutes. The tissue sections should not dry out during the pre-treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of FLEX IHC Microscope Slides (Code K8020) is recommended. After staining the sections must be dehydrated, cleared and mounted using permanent mounting medium.
Staining procedure including materials required but not supplied	The recommended visualization system is EnVision™ FLEX+, Mouse, High pH (Code K8012). The staining steps and incubation times are pre-programmed into the software of Dako Autostainer/Autostainer Plus instruments, using the following protocols: Template protocol: FLEXRTU2 (200 µL dispense volume) or FLEXRTU3 (300 µL dispense volume) Autoprogram: PR (without counterstaining) or Autoprogram: PRH (with counterstaining) The Auxiliary step should be set to "rinse buffer" in staining runs with ≤10 slides. For staining runs with >10 slides the Auxiliary step should be set to "none." This ascertains comparable wash times. All incubation steps should be performed at room temperature. For details, please refer to the Operator's Manual for the dedicated instrument. If the protocols are not available on the used Dako Autostainer instrument, please contact Dako Technical Services. Counterstaining in hematoxylin is recommended using EnVision™ FLEX Hematoxylin (Dako Autostainer/Autostainer Plus) (Code K8018). Non-aqueous, permanent mounting medium is recommended. Positive and negative controls should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include columnar epithelial cells, the basal squamous epithelial cells and the stromal cells and the cells/structures should display reaction patterns as described for this tissue in "Performance characteristics" in all positive specimens. The recommended negative control reagent is FLEX Negative Control, Mouse (Dako Autostainer/Autostainer Plus) (Code IS750).
Staining interpretation	The cellular staining pattern is nuclear. Cytoplasmic labeling, if observed, should be regarded as non-specific. A positive result is defined as nuclear staining in ≥1% of tumor cells at any stain intensity. This is consistent with ASCO/CAP's recommended cut-off of ≥1% positive tumor cells for positive assessment (11).
Product specific limitation	1. False-negative results could be caused by degradation of the antigen in the tissues over time. Specimens should be stained within 2 months of mounting of tissues on slides when stored at room temperature (12). 2. Lymphocytes and stromal cells can exhibit positive staining with the FLEX Monoclonal Mouse Anti-Human PR clone PgR 636 antibody and should not be interpreted or counted as PR staining tumor cells. 3. For optimal and reproducible results, the PR protein requires target retrieval when tissues are routinely fixed (neutral buffered formalin) and paraffin embedded. 4. Use of Dako Anti-PR, Clone PgR 636 on tissues with fixatives other than formalin has not been validated.
Performance characteristics	<u>Precision:</u> Serial sections from each of 12 different formalin-fixed paraffin-embedded blocks of breast carcinoma, representing a dynamic range of PR expression, were collected for testing. Testing was performed as follows: <u>Intra-run precision:</u> Following the standard EnVision™ FLEX+, High pH protocol, three sections from each tissue block were stained with Anti-PR, Clone PgR 636. Concurrently one section from each block was stained with a negative control reagent. <u>Inter-run precision:</u> Staining one section from each tissue block, the above procedure was repeated on five non-consecutive days. Concurrently, one section from each tissue block was stained with a negative control reagent. <u>Inter-instrument precision:</u> Staining a total of three sections from each tissue block, the above procedure was performed on three different Autostainer instruments by three different operators. Concurrently, one slide from each tissue block was stained with a negative control reagent. Precision experiments with Anti-PR, Clone PgR 636 yielded consistent results with intra-run, inter-run and inter-instrument testing. Consistent test conditions were maintained throughout the study and reagents were stored at 2-8 °C between test runs.

Normal tissues: Table 1 contains a summary of Anti-PR, Clone PgR 636 immunoreactivity on the recommended panel of normal tissues. All tissues were formalin-fixed and paraffin-embedded and stained with Anti-PR, Clone PgR 636 according to the instructions in the package insert. Cytoplasmic staining was observed with Anti-PR, Clone PgR 636 in several different tissue elements including epithelium, stroma, interstitial cells and inflammatory cells. While cytoplasmic staining was observed it is not considered diagnostic as per the intended use of this antibody.

Table 1: Summary of Anti-PR, Clone PgR 636 Normal Tissue Reactivity

Tissue Type (# tested)	Positive Tissue Elements	Tissue Type (# tested)	Positive Tissue Elements
Adrenal (3)	1/3 cells in glomerulosa region (50%), nuclear	Nerve, peripheral (3)	0/3
	1/3 cells in glomerulosa region (50%), nuclear	Ovary (3)	3/3 Stromal cells (50-70%), nuclear
Bone marrow (3)	0/3	Pancreas (3)	2/3 Islets of Langerhans (50-90%), nuclear
Breast (3)	2/3 Glandular epithelial cells (50-90%), nuclear	Parathyroid (3)	3/3 Glandular epithelial cells (1-10%), nuclear
Cerebellum (3)	0/3	Pituitary (3)	3/3 Pituitary glandular cells (1-40%), nuclear
Cerebrum (3)	1/3 Meningial cells (100%), nuclear	Prostate (3)	3/3 Stromal cells (30-80%), nuclear
Cervix (3)	3/3 Epithelial cells (50-90%), nuclear	Salivary gland (3)	3/3 Glandular epithelial cells (<1%-60%), nuclear
	3/3 Stroma, including inflammatory cells (50%), nuclear	Skin (3)	0/3
Colon (3)	1/3 Lymphoid/inflammatory cells (10%), nuclear	Small intestine (3)	3/3 Stromal and inflammatory cells (30-50%), nuclear
	1/3 Lymphoid/inflammatory cells (10%), nuclear	Spleen (3)	0/3
Esophagus (3)	1/3 Stromal cells (50%), nuclear	Stomach (3)	1/3 Interstitial cells (20%), nuclear
Kidney (3)	3/3 Interstitial cells (1-5%), nuclear	Testis (3)	3/3 Interstitial cells (5-80%), nuclear
Liver (3)	0/3	Thymus (3)	0/3
Lung (3)	2/3 Interstitial cells (1-10%), nuclear	Thyroid (3)	0/3
	2/3 Inflammatory cells (1-10%), nuclear	Tonsil (3)	0/3
Mesothelial cells (2)	0/2	Uterus (2)	2/2 Glandular epithelial cells (100%), nuclear
Muscle, cardiac (3)	3/3 Myocytes (30%), peri-nuclear		2/2 Myometrial stromal cells (100%), nuclear
Muscle, skeletal (3)	0/3		

Method comparison: Monoclonal Mouse Anti-Human PR, Clone PgR 636 testing was performed with EnVision™ FLEX+ and scored according to ASCO/CAP guidelines ($\geq 1\%$ cut-off) (11). Anti-PR (Clone 1294) testing was performed using Dako ER/PR pharmDx™ Kit and scored using the Allred scoring guideline described in the package insert. The method comparison data are presented in Table 2. Using these respective scoring guidelines, Monoclonal Mouse Anti-Human PR, Clone PgR 636 was concordant with the PR antibody component of Dako ER/PR pharmDx™ Kit, exhibiting values for overall, positive and negative agreement of 94.5%, 95.8% and 93.1%, respectively. Table 3 shows the comparison of both assays when scored using the ASCO/CAP guidelines.

Table 2: Agreement between Anti-PR, Clone PgR 636 (ASCO/CAP) and Anti-PR Component of ER/PR pharmDx™ Kit (Allred)

		Anti-PR Component of ER/PR pharmDx™ Kit		Total
		Positive	Negative	
Monoclonal Mouse Anti-Human PR, Clone PgR 636	Positive	115	8	123
	Negative	5	108	113
	Total	120	116	236

Positive Percent Agreement = 95.8% (95% CI: 91.1-96.8)

Negative Percent Agreement = 93.1% (95% CI: 90.5-96.7)

Overall Percent Agreement = 94.5% (95% CI: 90.8-96.8)

Table 3: Agreement between Anti-PR, Clone PgR 636 (ASCO/CAP) and Anti-PR Component of ER/PR pharmDx™ Kit (ASCO/CAP)

		Anti-PR Component of ER/PR pharmDx™ Kit		Total
		Positive	Negative	
Monoclonal Mouse Anti-Human PR, Clone PgR 636	Positive	115	8	123
	Negative	1	112	113
	Total	116	120	236

Positive Percent Agreement = 99.1% (95% CI: 93.0-98.0)

Negative Percent Agreement = 93.3% (95% CI: 92.8-98.0)

Overall Percent Agreement = 96.2% (95% CI: 93.0-98.0)

Inter Laboratory Reproducibility: Anti-PR, Clone PgR 636 reproducibility testing was performed in three testing laboratories over five non-consecutive days on 21 unique breast cancer specimens and scored according to ASCO/CAP guidelines ($\geq 1\%$ cut-off) for a total of 315 evaluations. Tables 4, 5, and 6 detail the site to site reproducibility of the assay. The average positive, average negative, and overall percent agreement calculations support the highly reproducible results of the PR (PgR 636) assay when used for the determination of PR status in a clinical setting.

Table 4: Site 1 vs. Site 2 Inter Laboratory Reproducibility of Anti-PR, Clone PgR 636

		Site 1		Total
Site 2		Positive	Negative	
	Positive	55	0	55
	Negative	4	46	50
	Total	59	46	105

Average Positive Percent Agreement = 96.5%

Average Negative Percent Agreement = 95.8%

Table 5: Site 1 vs. Site 3 Inter Laboratory Reproducibility of Anti-PR, Clone PgR 636

		Site 1		Total
Site 3		Positive	Negative	
	Positive	59	1	60
	Negative	0	45	45
	Total	59	46	105

Average Positive Percent Agreement = 99.2%

Average Negative Percent Agreement = 98.9%

Table 6: Site 2 vs. Site 3 Inter Laboratory Reproducibility of Anti-PR, Clone PgR 636

		Site 2		Total
Site 3		Positive	Negative	
	Positive	55	5	60
	Negative	0	45	45
	Total	55	50	105

Average Positive Percent Agreement = 95.7%

Average Negative Percent Agreement = 94.7%

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Bibliographie

Literaturnachweise

1. Henderson C. Breast cancer. In: Harrison's principles of internal medicine, 12th edition, Wilson JD, Braunwald E, Isselbacher KJ, Petersdorf RG, Martin JB, Fauci AS, Root RK (eds). McGraw-Hill, Inc. New York, 1991
2. Fuqua SAW. Estrogen and progesterone receptors and breast cancer. Diseases of the Breast, Harris et al, eds. Lippincott-Raven, 1996. p. 261.
3. McGuire WL, Clark GM. The prognostic role of progesterone receptors in human breast cancer. Sem Oncol 1983;10:2.
4. Clark GM, McGuire WL, Hubay CA, Pearson OH, Marshall JS. Progesterone receptors as prognostic factor in stage II breast cancer. N Engl J Med 1983;39:1343.
5. Ravdin PM, Green S, Dorr TM, McGuire WL, Fabian C, Puch RP, et al. Prognostic significance of progesterone receptor levels in estrogen receptor-positive patients with metastatic breast cancer treated with tamoxifen: Results of a prospective southwest oncology group study. JCO 1992;10:1284.
6. Chevallier B, Heintzmann F, Mosseri V, Dauce JP, Bastit P, Graic Y, Brunelle P, Basuyay JP, Comoz M, Asselain B. Prognostic value of estrogen and progesterone receptors in operable breast cancer: Results of a univariate and multivariate analysis. Cancer 1988; 62:2517
7. Page DL, Jensen RA, Simpson JF. Routinely available indicators of prognosis in breast cancer. Breast Can Res Treat 1998;51:195.
8. Allred DC, Harvey JM, Berardo M, Clark GM. Prognostic and predictive factors in breast cancer by immunohistochemical analysis. Mod Pathol 1998;11:155.
9. Fitzgibbons PL, Page DL, Weaver D, Thor AD, Allred DC, Clark GM, et al. Prognostic factors in breast cancer. College of American Pathologists consensus statement. Arch Pathol Lab Med 2000;124:966.
10. Press M, Spaulding B, Groshen S, Kaminsky D, Hagerty M, Sherman L, et al. Comparison of different antibodies for detection of progesterone receptor in breast cancer. Steroids 2002; 67:799.
11. Hammond MEH, Hayes DF, Dowsett M, Allred C, Hagerty KL, Badve S, et al. American Society of Clinical Oncology/College of American Pathologists Guideline Recommendations for Immunohistochemical Testing of Estrogen and Progesterone Receptors in Breast Cancer. Arch Pathol Lab Med 2010;134:907–22.
12. Clinical and Laboratory Standards Institute (formerly NCCLS). Quality assurance for immunocytochemistry; Approved guideline. CLSI document MM4-A (1-56238-396-5)- CLSI, 940 West Valley Road, Suite 1400. Wayne, PA 19087-1898 USA 1999.

 Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	 Dako North America, Inc. 6392 Via Real Carpinteria, California 93013 USA Tel 805 566 6655 Fax 805 566 6688 Technical Support 800 424 0021 Customer Service 800 235 5763	 Dako Denmark A/S Produktionsvej 42 DK-2600 Glostrup Denmark Tel +45 4485 9500 Fax +45 4485 9595 www.dako.com
 Manufacturer Fabricant Hersteller	 Batch code Code du lot Chargenbezeichnung	 Contains sufficient for <N> tests Contenu suffisant pour <n> tests Inhalt ausreichend für „n“ Ansätze		
 Use by Utiliser jusque Verwendbar bis	 Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	 Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU		

PT0020/ Rev C

Edition / Édition / Ausgabe 08/14